
Timberline Series

Electronic Precision Balances/Elektronische
Präzisionswaagen/Balances électroniques
de précision/Balanzas electrónicas de precisión

Operating Instructions/Betriebsanleitung/
Mode d'emploi/Instrucciones de funcionamiento

English – page 3

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Contents

Warnings and Safety Precautions	4
Getting Started	5
Operation	10
Basic Weighing Function	10
Calibration/Adjustment	11
Application Programs	12
Net-Total Formulation/Second Tare Memory	12
Counting	13
Weighing in Percent	14
Weigh Averaging	15
Toggling between Weight Units	16
Configuring the Balance	17
Setting the Parameters (Menu Codes)	17
Balance Operating Menu (Overview)	18
ISO/GLP-compliant Printout	20
Data Interface	22
Error Codes	23
Care and Maintenance	24
Instructions for Recycling	26
Overview	27
Specifications	27

Warnings and Safety Precautions

Safety

- To prevent damage to the equipment, please read these operating instructions carefully before using your balance.
- ⚠ Do not use this balance in a hazardous area/location.
- ⚠ Make absolutely sure to unplug the balance from AC power before you connect or disconnect a peripheral device.
- ⚠ Exposure to excessive electromagnetic disturbance can cause the readout value to change. Once the disturbance has ceased, the instrument can be used again in accordance with its intended use.

Setting up the Balance

- ⚠ Warning when using pre-wired RS-232 connecting cables: The pin assignments in RS-232 cables purchased from other manufacturers may be incompatible with DENVER balances. Be sure to check the pin assignment against the chart on page 21 before connecting the cable, and disconnect any lines that do not match.
- Connect only DENVER accessories and options, as these are optimally designed for use with your DENVER balance. Do not try to solve any problems on your own. The operator shall be responsible for any modifications to DENVER equipment and for any connections of cables or equipment not supplied by DENVER and must check and, if necessary, correct these modifications and connections. On request, DENVER will provide information on the minimum operating specifications (in accordance with the standards for defined immunity to interference).
- Do not open the balance housing. If the seal is broken, this will result in forfeiture of all claims under the manufacturer's warranty.

Getting Started

Storage and Shipping Conditions

- Do not expose the balance to extreme temperatures, blows, shocks, vibration or moisture.

Unpacking the Balance

- After unpacking the balance, check it immediately for any visible damage.
- If you see any sign of damage, proceed as directed in the chapter entitled “Care and Maintenance,” under the section on “Safety Inspection.”
- Save the box and all parts of the packaging until you have successfully installed your balance in case you need to return it. Before packing your balance, unplug all connected cables to prevent damage.

Equipment Supplied

- Balance
- Weighing pan
- Pan support (only on models with a round weighing pan)
- AC adapter, plug type

Additionally supplied with TP-214, TP-114, TP-64 and TP-323DS, TP-153DS models:

- Shield ring (not for models TP-...DS)
- Shield plate
- Dust cover

Additionally supplied with TP-323 and TP-153 models:

- Glass draft shield with cover


Installation Instructions

When choosing a location to set up your balance, observe the following:

- Avoid placing the balance in close proximity to a heater or otherwise exposing the balance to heat or direct sunlight
- Protect the balance from drafts that come from open windows or doors
- Avoid exposing the balance to extreme vibrations during weighing
- Do not expose the balance to extreme moisture over long periods

Conditioning the Balance


Moisture in the air can condense on the surfaces of a cold balance whenever it is brought into a substantially warmer place. If you transfer the balance to a warmer area, make sure to condition it for about 2 hours at room temperature, leaving it unplugged from AC power.

Setting up the Balance


Balances with an Analytical Draft Shield

- Place the components listed below inside the chamber in the order given:
 - Shield plate
 - Shield ring (not for TP-...DS models)
 - Pan support
 - Weighing pan

Balances with a Glass Draft Shield

- Place the components listed below inside the chamber in the order given:
 - Draft shield base – place it on the balance so that the edge for fitting the glass draft shield faces upwards and turn it until it is firmly in place
 - Pan support
 - Weighing pan
 - Glass draft shield
 - Draft shield cover – place it on the balance so that the edge faces downwards


Balances with a Round Weighing Pan

- Place the components listed below inside the chamber in the order given:
 - Pan support
 - Weighing pan

Balances with a Rectangular Weighing Pan

- Place the weighing pan on the balance


Connecting the Balance to AC Power/ Safety Precautions

Use only original Denver AC adapters


- Insert the right-angle plug into the jack
- The AC adapter rated to Class 2 can be plugged into any wall outlet without requiring any additional safety precautions

The ground is connected to the balance housing, which can be additionally grounded for operation.

Using a Non-Rechargeable/Rechargeable Battery: for models TP-6101, TP-3101, TP-12

- A non-rechargeable or rechargeable battery is not included with the equipment supplied
- ⚠ Use only a commercially available non-rechargeable or rechargeable 9-volt battery
- ⚠ When using a rechargeable battery, always use an external charger to recharge the battery
- Lay the balance/scale on its side
- Open the battery compartment:
Lift the compartment cover

- Insert the 9-volt (rechargeable) battery in the compartment
- Make sure to connect the positive and negative poles correctly
- ⚠ All used batteries are classified as waste that requires special handling (not "household" waste). Dispose of rechargeable batteries in accordance with the applicable special waste disposal regulations.
- Close the battery compartment:
Press down on the cover until it clicks into place

Federal Communications Commission (FCC) Statement

NOTE: This equipment has been tested and found to comply with the limits pursuant to part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference to radio communications.

For information on the specific limits and class of this equipment, please refer to the Declaration of Conformity. Depending on the particular class, you are either required or requested to correct the interference.


If you have a Class A digital device, you need to comply with the FCC statement as follows:
“Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.”

If you have a Class B digital device, please read and follow the FCC information given below:

“However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.”

Before you operate this equipment, check which FCC class (Class A or Class B) it has according to the Declaration of Conformity included. Be sure to observe the information of this Declaration.


Leveling the Balance

Level the balance any time you set it up in a new location. Use only the 2 front feet of the balance for leveling.

TP-3102, TP-1502, TP-6101, TP-3101, TP-12:

- Turn the 2 rear feet until they are in position
 - Turn the 2 front feet as shown here in the illustration until the air bubble is centered in the level indicator
- > In most cases, this will require several adjustment steps

Anti-theft Locking Device

To protect against theft, use the mounting lug on the rear panel of the balance.

- Secure the balance at the place of installation, for example with a chain or a lock

Operation


Basic Weighing Function

Preparation

- Turn on the balance:
Press [ON/OFF]
- To change configurations: See the chapter entitled "Configuring the Balance"
- To tare the balance: Press [TARE]

Additional Functions:

- To turn off the balance:
Press [ON/OFF]

Warmup Time

- To ensure accurate results, the balance must warm up for 60 minutes before operation. Only after this time will the balance have reached the required operating temperature.

Example

Basic weighing

Step	Key (or instruction)	Display/Printout
1. Turn on the balance	[ON/OFF]	
Self-test is performed		
2. Place container on balance (here: 52.0 g)		+ 52.0 g
3. Tare the balance	[TARE]	+ 0.0 g
4. Place sample in container on balance (here: 150.2 g)		+ 150.2 g

Calibration and Adjustment

Calibration is the determination of any difference between the weight readout and the true weight (mass) of a sample. Adjustment is the correction of this difference, or its reduction to an allowable level within maximum permissible error limits.

In the Timerline Series balance, both functions are performed at one touch of the CAL key. This is why we use the term “calibrate/adjust” in the example given below.

Available Features

Calibration/adjustment can only be performed when

- there is no load on the balance,
- the balance is tared,
- the internal signal is stable.



If these conditions are not met, an error message is displayed. Otherwise, the weight required for calibration/adjustment is displayed.

You can use any of the following weight units to calibrate/adjust:
g, kg*, lb (menu code 1.4.x)

You can block calibration/adjustment of the balance:

- Select menu code 1.5.3

External calibration/adjustment

Step	Key (or instruction)	Display/Printout
1. Zero the balance	[TARE]	0.0 g
2. Begin calibration/adjustment Calibration weight is displayed without weight unit (here: 1000 g)	[CAL]	+ 1000.0
3. Place the indicated calibration weight on the balance		1000.0
After calibration/adjustment, the calibration weight is displayed with wt. unit		+ 1000.0 g
4. Remove the calibration weight		0.0 g

* = not available for models with a readability of 0.1 mg

Application Programs




Net-Total Formulation/Second Tare Memory

With this application program you can weigh in components for formulation of a mixture.

Preparation

Configure the Net-Total Formulation/Second Tare Memory application in the operating menu: See "Configuring the Balance." Menu code: 2.13

Example

Step	Key (or instruction)	Display/Printout
1. Place an empty container on the balance		+ 65.0 g
2. Tare the balance	[TARE]	+ 0.0 g
3. Add the first component		+ 120.5 g
4. Store the first component weight. If the print format is set to include data ID codes, the following is printed	[FUNCTION]	0.0 g _{NET}
		N1 + 120.5 g
5. Add the next component		+ 70.5 g
6. Store the 2nd component weight.	[FUNCTION]	0.0 g _{NET}
7. Add further components, if desired	As described for steps 5 and 6	
8. Display total weight and fill to desired final weight	[CF]	+ 191.0 g

Counting

Purpose

With the Counting program you can determine the number of parts that each have approximately equal weight.

○ Reference sample quantity:

- Code 3.3.1 5 pcs
- Code 3.3.2 10 pcs (factory setting)
- Code 3.3.3 20 pcs
- Code 3.3.4 50 pcs
- Code 3.3.5 100 pcs

Preparation

- Configure the Counting application in the operating menu:
See "Configuring the Balance"
Menu code: 2.1.4

- Storage parameter (display accuracy for counting)
Code 3.4.1 Standard resolution (factory setting)
Code 3.4.2 With 10 times higher resolution than standard




See also "Configuring the Balance"

Example

Determine an unknown piece count; weigh the preset reference sample quantity

Menu: Application program: Counting (menu code 2.1.4)

Reference sample quantity: 20 pcs (menu code 3.3.3)

Step	Key (or instruction)	Display/Printout
1. Tare the balance	[TARE]	0.0 g
2. Display the reference sample quantity (here: 20 pcs)	[FUNCTION] >2 sec.	rEF 20 (briefly)
3. Place the reference sample quantity (20 pcs) on the balance (here: 66 g)		+ 66.0 g
4. Start the application; if the print format is set to include data ID codes, the following piece weight is printed	[FUNCTION]	+ 20 pcs wRef + 3.300 g
5. Weigh uncounted parts (here: 174 pcs)		+ 174 pcs
6. Display weight	[FUNCTION]	+ 574.2 g
7. Display quantity	[FUNCTION]	+ 174 pcs
8. Unload the balance		0 pcs
9. Delete the reference value	[CF]	
10. Repeat the procedure starting from step 5, if desired.		

Weighing in Percent

Purpose

This application program allows you to obtain weight readouts in percent which are in proportion to a reference weight.

- Reference percentage:
 - Code 3.3.1 5 %
 - Code 3.3.2 10 % (factory setting)
 - Code 3.3.3 20 %
 - Code 3.3.4 50 %
 - Code 3.3.5 100 %

Preparation




- Configure the Weighing in Percent application in the operating menu: See "Configuring the Balance."
Menu code: 2.1.5
- Storage parameter (display accuracy for counting)
 - Code 3.4.1 Standard resolution: With stability (factory setting)
 - Code 3.4.2 With 10 times higher stability than standard
 See also "Configuring the Balance"

Example

Determine an unknown percentage: store the weight on the balance as a reference percentage

Menu: Application program: Weighing in Percent (menu code 2.1.5)

Menu: Reference percentage: 100 % (menu code 3.3.5)

Step	Key (or instruction)	Display/Printout
1. Tare the balance	[TARE]	0.0 g
2. Display the reference percentage:	[FUNCTION] >2 sec.	rEF 100
3. Place the reference weight for 100 % on the balance (here: 222.5 g)		+ 222.5 g
4. Start the application; if the print format is set to include data ID codes the following is printed	[FUNCTION]	+ 100.00 % Wxx% + 222.500 g
5. Place an unknown weight on the balance (here: 322.5 g)		+ 144.94 %
6. Display weight	[FUNCTION]	+ 322.5 g
7. Display percentage	[FUNCTION]	+ 144.94 %
8. Unload the balance		0.00 %
9. Delete the reference percentage	[CF]	
10. Repeat the procedure starting from step 5, if desired.		

Weigh Averaging

Purpose

Use this program to determine weights under unstable ambient conditions. In this program, the balance calculates the weight as the average value from a defined number of individual weighing operations, called "subweighs" in this manual.

- Number of subweighs for weigh averaging:

3.3.1	5 subweighs
3.3.2	10 subweighs (factory setting)
3.3.3	20 subweighs
3.3.4	50 subweighs
3.3.5	100 subweighs

Preparation



- Configure the Weigh Averaging application in the operating menu: See "Configuring the Balance."
Menu code: 2.1.12

See also "Configuring the Balance"

Example

Determine the weight of a sample in extremely unstable ambient conditions by calculating the average of 10 subweighing operations.

Menu: Application program: Weigh Averaging (menu code 2.1.12)

Step	Key (or instruction)	Display/Printout
1. Tare the balance	[TARE]	0.0 g
2. Display the number of subweighs (here: 10)	[FUNCTION] >2 sec.	rEF 10 (briefly)
3. Place sample on the balance (weight readout fluctuates)		8888
4. Start measurement	[FUNCTION]	8888 10 9 8 ⋮ 1
After 10 subweighs		+ 275.5 g Δ
If the print format is set to include data ID codes, the following is printed		Res + 275.5 g
5. Unload the balance		+ 275.5 g Δ (stable display)
6. Delete the result	[CF]	
7. Repeat the procedure starting from step 3, if desired.		

Toggleing between Weight Units

With this application program you can toggle the display of a weight value back and forth between two weight units.

Configure the “Toggle Weight Units” application in the operating menu:
See “Configuring the Balance.” Menu code 2.1.2

Menu code		Unit	Conversion factor	Abbr. on printout
1.7.1	3.1.1	Grams	1.0000000000	o
1.7.2 o	3.1.2 o	Grams	1	g
1.7.3	3.1.3	Kilograms ¹⁾	0.0010000000	kg
1.7.4	3.1.4	Carats	5	ct
1.7.5	3.1.5	Pounds	0.00220462260	lb
1.7.6	3.1.6	Ounces	0.03527396200	oz
1.7.7	3.1.7	Troy ounces	0.03215074700	ozt
1.7.8	3.1.8	Hong Kong taels	0.02671725000	tlh
1.7.9	3.1.9	Singapore taels	0.02645544638	tls
1.7.10	3.1.10	Taiwanese taels	0.02666666000	tlt
1.7.11	3.1.11	Grains	15.43235835000	GN
1.7.12	3.1.12	Pennyweights	0.64301493100	dwt
1.7.13	3.1.13	Milligrams	1000	mg
1.7.14	3.1.14	Parts per pound	1.12876677120	/lb
1.7.15	3.1.15	Chinese taels	0.02645547175	tlc
1.7.16	3.1.16	Mommes	0.26670000000	mom
1.7.17	3.1.17	Austrian carats	5	K
1.7.18	3.1.18	Tola	0.08573333810	tol
1.7.19	3.1.19	Baht	0.06578947437	bat
1.7.20	3.1.20	Mesghal	0.21700000000	MS

o = Factory setting

¹⁾ = not for models with a readability of ≤ 0.1 mg

Function

- To toggle the display between the 1st and 2nd weight units:
Press the [FUNCTION] key

Configuring the Balance

Setting the Parameters (Menu Codes)

You can configure your balance to meet individual requirements by selecting from the parameters available in the menu.


Example: Adapt the balance to unstable ambient conditions
Menu code 1. 1. 4

Step	Key (or instruction)	Display
1. Turn off the balance	[ON/OFF]	
2. Turn the balance back on; while all segments are displayed	[ON/OFF]	1.
<input type="radio"/> To navigate within a menu level; the last menu option is followed by the first option	[TARE] repeatedly	2. : 9. 1.
3. Select the 2nd menu level	[PRINT]	1. 1.
4. Select the 3rd menu level	[PRINT]	1. 1. 2 o
5. In Menu Level 3: Select the desired option	[TARE] repeatedly	1. 1. 4
6. Confirm new setting; the "o" indicates the currently set option	[PRINT] >2 sec.	1. 1. 4 o
<input type="radio"/> Select the next menu level (here: move from the 3rd to the 1st level)	[PRINT]	1.
<input type="radio"/> Set other menu codes, if desired	[PRINT], [TARE]	
7. Store parameter settings and exit operating menu or	[TARE] >2 sec.	
<input type="radio"/> Exit operating menu without storing changes	[ON/OFF]	
> Restart the application		0.0 g


Balance Operating Menu (Overview)

o Factory setting

✓ User setting

¹⁾ = not for models with a readability of 0.1 mg

ISO/GLP-compliant Printout

Features

You can have the parameters pertaining to the ambient weighing conditions printed before (GLP header) and after (GLP footer) the values of a weighing series.

These parameters include:

GLP header:

- Date
- Time at beginning of measurement
- Balance manufacturer
- Balance model
- Balance serial number
- Software version number
- Identification number of the current sampling operation

GLP footer:

- Date
- Time at end of measurement
- Field for operator signature

△ The record can only be output to a Denver data printer 901042-1

Settings

- Set the following menu codes (see “Configuring the Balance”):
 - GLP-compliant record: menu code 5 5 2
 - ISO/GLP-compliant record after calibration/adjustment only: menu code 7 2 2 or ISO/GLP-compliant record always on: menu code 7 2 3
 - Line format for printout: With data ID codes – 22 characters: menu code 7 1 2
 - Printout date/time:
 - 24-h format: menu code 7 3 1
 - 12-h format: menu code 7 3 2
- △ No ISO/GLP-compliant record is output if any of the following settings are configured: menu codes 5 1 3, 5 1 4 (automatic printout) and 7 1 1

Function Keys

Press [PRINT] to output header and first measured value.

- > Header is output the first time [PRINT] is pressed

To output header and reference data automatically with an application program active: Press [FUNCTION]

End application program:
End application program and output GLP footer: Press [CF]

The ISO/GLP-compliant record can contain the following lines:

-----			Dotted line
17-Jan-2007		10:15	Date/time (beginning of measurement)
Denver	Inst		Balance manufacturer
Mod.		TP-3102	Balance model
Ser. no.		10105355	Balance serial number
Ver. no.		00-19-41	Software version
ID			ID
-----			Dotted line
L ID			Measurement series no.
wRef	+	21.14 g	Counting: Reference weight
Qnt	+	235 pcs	Counting result
Qnt	+	567 pcs	Counting result
-----			Dotted line
17-Jan-2007		10:20	Date/time (end of measurement)
Name:			Field for operator signature
-----			Blank line
-----			Dotted line

ISO/GLP-compliant printout for external calibration/adjustment


-----			Dotted line
17-Jan-2007		10:30	Date/time (beginning of measurement)
Denver	Inst		Balance manufacturer
Mod.		TP-3102	Balance model
Ser. no.		10105355	Balance serial number
Ver. no.		00-19-41	Software version
ID			ID
-----			Dotted line
Cal. Extern			Calibration/adjustment mode
Set	+	2000.0 g	Calibration weight
-----			Dotted line
17-Jan-2007		10:32	Date/time (end of measurement)
Name:			Field for operator signature
-----			Blank line
-----			Dotted line

Data Interface

Purpose

Your balance comes equipped with an interface port for connection to a computer or other peripheral device. You can use an on-line computer to change, start and/or monitor the functions of the balance and the application programs.

Female interface connector


Pin Assignment Chart, 25-pin female interface connector, RS-232:

- Pin 1: Shield
- Pin 2: Data output (TxD)
- Pin 3: Data input (RxD)
- Pin 4: Internal ground (GND)
- Pin 5: Clear to send (CTS)
- Pin 6: Not connected
- Pin 7: Internal ground (GND)
- Pin 8: Internal ground (GND)
- Pin 9: Not connected
- Pin 10: Not connected
- Pin 11: Charging voltage for rechargeable battery pack +12 ... +20 V (I_{out} 25mA)
- Pin 12: Reset _ Out *)
- Pin 13: +5 V output
- Pin 14: Internal ground (GND)
- Pin 15: Universal remote switch
- Pin 16: Not connected
- Pin 17: Not connected
- Pin 18: Not connected
- Pin 19: Not connected
- Pin 20: Data terminal ready (DTR)
- Pin 21: Internal ground (GND)
- Pin 22: Not connected
- Pin 23: Not connected
- Pin 24: Not connected
- Pin 25: +5 V output

Preparation

You can set these parameters for other devices in the Setup menu (see the chapter entitled "Configuring the Balance"). You will also find a detailed description of the available data interface commands in the file "Data Interface Descriptions for TP Models", which you can order directly from Denver Instrument.

The many and versatile properties of these balances can be fully utilized for printing out records of the results when you connect your balance to a DENVER data printer. The recording capability for printouts makes it easy for you to work in compliance with ISO/GLP.

*) = Hardware restart

Error Codes

Error codes are shown on the main display for 2 seconds. The program then returns automatically to the previous mode (e.g., weighing).

Display	Cause	Solution
No segments appear on the display	No AC power is available The AC adapter is not plugged in	Check the AC power supply Plug in the AC adapter
	Battery or rechargeable battery pack is discharged	Replace the battery or recharge the battery pack using an external charger
H	The load exceeds the balance capacity	Unload the balance
L and E 54	The weighing pan is not in place Something is touching the weighing pan	Place the weighing pan on the balance Move that object that is touching the weighing pan
E 02	Calibration parameter not met, e.g.: – balance not zeroed – balance is loaded	Unload the balance Press [TARE] to tare the balance Calibrate only when zero is displayed
E 09	When gross value \leq zero; no tare	Tare the balance
E 10	The [TARE] key is blocked when there is data in the second tare memory (net-total). Only 1 tare function can be used at a time	Press [CF] to clear the tare memory and release the tare key
E 11	Value input is not allowed for second tare memory	Press [TARE]
E 22	Weight is too light or there is no sample on the balance	Increase the weight on the balance
E 30	Interface port for printer output is blocked	Contact your local DENVER Service Center
Max. weighing capacity is less than indicated under “Specifications”	The balance was turned on without the weighing pan in place	Place the weighing pan on the balance and press [ON/OFF] to turn the balance back on
The weight readout is obviously wrong	The balance has not been calibrated/adjusted The balance was not tared before weighing	Calibrate/adjust the balance Tare the balance

If any other errors occur, contact your local DENVER Service Center!


Care and Maintenance

Service

Regular servicing by a trained DENVER dealer will extend the service life of your balance and ensure its continued weighing accuracy. DENVER can offer you service contracts, with your choice of regular maintenance intervals. The optimum maintenance interval depends on the operating conditions at the place of installation and on the individual tolerance requirements.

Repairs

Repair work must be performed by trained service technicians. Any attempt by untrained persons to perform repairs may lead to hazards for the user.

Cleaning

- Unplug the AC adapter from the wall outlet (mains supply). If you have an interface cable connected to the balance port, unplug it from the port,
- Clean the balance using a piece of cloth which has been wet with a mild detergent (soap)
- After cleaning, wipe down the balance with a soft, dry cloth.

Removing and Cleaning the Weighing Pan:

- Lift up and remove the weighing pan together with the pan support by gripping them from under the shield ring. Make sure that you do not damage the weighing system in doing so.
- ⚠ Make sure that no liquid enters the balance housing.
- ⚠ Do not use any aggressive cleaning agents (solvents or similar agents).

Cleaning Stainless Steel Surfaces

Clean all stainless steel parts regularly. Remove the stainless steel weighing pan and thoroughly clean it separately. Use a damp cloth or sponge to clean any stainless steel parts on the balance. Only use commercially available household cleaning agent that is suitable for use on stainless steel. Clean stainless steel surfaces by wiping them down. Then rinse thoroughly, making sure to remove all residues. Afterwards, allow the balance to dry. If desired, you can apply oil to the cleaned surfaces as additional protection. Solvents are permitted for use only on stainless steel parts.

Safety Inspection

If there is any indication that safe operation of the balance with the AC adapter is no longer warranted:

- Turn off the power and disconnect the equipment from AC power immediately
- > Lock the equipment in a secure place to ensure that it cannot be used for the time being

In such cases, notify your DENVER service center. Maintenance and repair work may only be performed by service technicians who are authorized by DENVER and who

- have access to the required maintenance manuals
- have attended the relevant service training courses

Instructions for Recycling

Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste.

The EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it. In Germany and many other countries, Denver Instrument takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center:

Denver Instrument GmbH
Robert-Bosch-Breite 10
37079 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Denver affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Denver Instrument, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet/manual or visit our Internet website (www.denverinstrument.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

Overview

Specifications

Model		TP-214	TP-114	TP-64
Weighing capacity	g	210	110	60
Readability	mg	0.1	0.1	0.1
Tare range (subtractive)	g	210	110	60
Repeatability	±mg	0.1	0.1	0.1
Linearity	±mg	0.2	0.2	0.2
Allowable ambient operating temperature		+10 to +30°C (50°F to 86°F)		
Sensitivity drift within +10 to +30°C	±%/K	3 · 10 ⁻⁶		
Response time (average)	s	3		
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels		
Display update (depends on the filter level selected)	s	0.2–0.4		
External calibration weight (of at least accuracy class...)	g lb	200 (E2) 0.4	100 (E2) 0.2	50 (E2) 0.1
Net weight, approx.	kg/lb	3.0/6.6		
Pan size	mm inches	80 Ø 3.3 Ø		
Weighing chamber height	mm inches	200 7.9		
Dimensions (W×D×H)	mm inches	189×251×299 7.4×9.9×11.8		
AC power source/power requirements		AC adapter, 230 V or 115 V, +15% to -20%		
Frequency		48–60 Hz		
Power requirements, direct current	V	10 to 20		
Power consumption (average)	W	1		
Hours of operation with the rechargeable battery pack (see “Accessories”)	h	20		

Model		TP-323, TP-323DS	TP-153, TP-153DS	TP-3102	TP-1502
Weighing capacity	g	320	150	3,100	1,500
Readability	g	0.001	0.001	0.01	0.01
Tare range (subtractive)	g	320	150	3,100	1,500
Repeatability	≤±g	0.001	0.0015	0.01	0.015
Linearity	≤±g	0.002	0.003	0.02	0.03
Allowable ambient operating temperature		+10 to +30°C (50° to 86°F)			
Sensitivity drift within +10 to +30°C ≤±/K		4 · 10 ⁻⁶	4 · 10 ⁻⁶	3 · 10 ⁻⁶	4 · 10 ⁻⁶
Response time (average)	s	3	2.5	3	2,5
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels			
Display update (depends on the filter level selected)	s	0.2–0.8			
External calibration weight (of least accuracy class...)	g lb	200 (E2) 0.4	100 (F2) 0.2	2,000 (E2) 4	1,000 (F2) 2
Net weight, approx.	kg	2.2	2.2	2.2	1.6
Pan size	mm	100 ∅	100 ∅	174×143	174×143
Dimensions (W×D×H)	mm	189× 251×120 189× 251×299	189× 251×120 189× 251×299	189× 251×70	189× 251×70
AC power source/ power requirements		AC adapter 230 V or 115 V, +15% to -20%			
Frequency		48–60 Hz			
AC power source/direct current	V	10–20			
Power consumption (average)	W	0.75			
Hours of operation with the rechargeable battery pack (see "Accessories")	h	20	20	20	25

Model		TP-6101	TP-3101	TP-12
Weighing capacity	g	6,100	3,100	12,000
Readability	g	0.1	0.1	1
Tare range (subtractive)	g	6,100	3,100	12,000
Repeatability	≤±g	0.1	0.1	1
Linearity	≤±g	0.2	0.2	2
Operating temperature range		+10 to +30°C (50°F to 86°F)		
Sensitivity drift within +10 to +30 °C	≤±/K	5 · 10 ⁻⁶	10 · 10 ⁻⁶	25 · 10 ⁻⁶
Response time (average)	s	2	2	1.5
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels		
Display update (depends on the filter level selected)	s	0.2–0.8		
External calibration weight (of at least accuracy class)	kg lb	5 (F2) 10	1 (F2) 2	5 (M1) 10
Net weight, approx.	kg	1.6		
Pan size	mm	174×143		
Dimensions (B×T×H)	mm	189×251×70		
AC power source/, power requirements		AC adapter 230 V or 115 V, +15% to –20%		
Frequency		48–60 Hz		
AC power source/direct current	V	10–20		
Power consumption (average)	W	0.75		
Approx. hours of operation with – 9-volt battery, approx. – fully charged	h	11		
rechargeable battery	h	2.5		
Hours of operation with the YRB08Z rechargeable battery pack	h	25		

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Copyright by Denver Instrument.

All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior written permission of Denver Instrument.

The status of the information, specifications and illustrations in this manual is indicated by the date given below. Denver Instrument reserves the right to make changes to the technology, features, specifications, and design of the equipment without notice.

Status: October 2007, Denver Instrument.

Europe, Asia and Australia:

**Denver Instrument GmbH
Robert-Bosch-Breite 10
37079 Goettingen, Germany
Tel: +49-551-20977-30
Fax: +49-551-20977-39**

Internet:

www.denverinstrument.com

U.K. and Ireland:

**Denver Instrument Company
Denver House, Sovereign Way
Trafalgar Business Park
Downham Market
Norfolk PE38 9SW England
Tel: +44-136-63862-42
Fax: +44-136-63862-04**

North and South America:

**Denver Instrument Company
1401 17th Street, Suite 750
Denver, Colorado 80202
1-800-321-1135
Tel: +1-303-431-7255
Fax: +1-303-423-4831**