



## *Alkalinity of Bleach (NaClO Solution)*

### **Introduction**

A bleach solution contains mainly NaClO plus basic products such as NaOH and Na<sub>2</sub>CO<sub>3</sub>. The concentration of basic products is generally expressed as NaOH concentration, which is around 8 g/l NaOH, i.e. 0.2N concentration. The alkalinity determination uses an acid/base titration.

### **Principle**

The OH-content is simply determined by an acid/base titration using a 0.1 eq/l strong acid as titrant.

Depending on the bleach solution, the titration curve generally shows 2 inflections. The most common method involves a titration of all the basic functions by a predetermined end point titration at pH 4.00. Before this titration, it is necessary to reduce the ClO<sup>-</sup> ions present in the solution.

### **Electrodes and reagents**

pHC2011-8 Combined pH Electrode (part no. E16M317)  
 H<sub>2</sub>O<sub>2</sub> 30% (volume)

Hydrochloric acid 0.1 eq/l (or 0.1 mol/l)

Slowly add 8.3 ml of concentrated hydrochloric acid to 500 ml of distilled water and dilute to exactly 1000 ml. Calibrate the titrant versus Na<sub>2</sub>B<sub>4</sub>O<sub>7</sub> · 10 H<sub>2</sub>O (sodium borate) as standard.

(See separate application note)

Distilled water

IUPAC Series pH standards  
 pH 4.005 (part no. S11M002) and pH 10.012 (part no. S11M007)

### **End Point titration settings**

|                       |               |
|-----------------------|---------------|
| Burette volume:       | 25 ml         |
| Stirring speed:       | 400 rpm       |
| Working mode:         | pH            |
| Number of end points: | 1             |
| End point:            | 4.00 pH       |
| Stirring delay:       | 45 seconds    |
| Minimum speed:        | 0.2 ml/min    |
| Maximum speed:        | 8.0 ml/min    |
| Proportional band:    | 3.00 pH       |
| End point delay:      | 5 seconds     |
| Sample unit:          | ml            |
| Sample amount:        | 5             |
| Titration:            | Decreasing pH |
| Result:               | g/l           |

### **Procedure**

Calibrate the electrode using the two IUPAC standards

Pipette 5 ml of sample

Always dilute the sample with the same volume of distilled water (max. 50 ml)

Slowly add 5 ml of H<sub>2</sub>O<sub>2</sub> for 5 ml of sample

Dip electrode and delivery tip in the solution

Start method by pressing the RUN key

### **Results**

Expressed as NaOH content (MW = 40 g/mol) in g/l

As 1 molecule of titrant reacts with 1 molecule of NaOH

$$R(\text{NaOH}) = V(\text{titr}) * C(\text{titr}) * 40 / V(\text{smp})$$

-V(titr) = total volume of titrant to reach the end point in ml

-V(smp) = sample amount

-C(titr) = exact concentration of the titrant in eq/l

### **For a result in g/l**

Enter

The sample amount in the SAMPLE screen

The titrant concentration in the TITRANT screen

1 Titrant and 1 Sample in the COEFFICIENTS display

40 for NaOH molecular weight

The Titration Manager gives a result according to the above formula.

### **For a result as a %**

As the Titration Manager cannot give a result in % if the sample unit is a volumetric unit, you can use the equation feature:

Equation number: 1  
 Equation result: % NaOH  
 Equation formula: R1 / 10  
 R1 is the titration result calculated in g/l.

### **5 determinations on a commercial concentrated bleach**

Mean (as NaOH): 7.25 g/l  
 Standard deviation: 0.07 g/l  
 Rel. standard deviation: 1%

### **Working Range**

---

Results are expressed in g/l of NaOH (MW = 40 g/mol)

1 ml of HCl 0.1 eq/l represents 4 mg of NaOH or, with a 5 ml sample volume a NaOH content corresponding to 0.8 g/l

The working range can be calculated as the following formula:

$$\text{Result (in g/l)} = V(\text{titr in ml}) * 0.8 * 1000 / V(\text{smp})$$

Using the conditions given in this application note (5 ml sample and a 25 ml burette for titrant), it is possible to obtain results between 7 g/l (for 35% capacity of the burette) and 20 g/l (total capacity of the burette) with the best possible accuracy and reproducibility.

### **Notes**

---

Addition of H<sub>2</sub>O<sub>2</sub> is necessary to reduce the ClO<sup>-</sup> ion to Cl<sup>-</sup> before running the alkalinity titration.