


XL-SERIES TOPLOADING BALANCES

FEATURES

- RS-232C bi-directional interface
- Multiple weighing units - grams, ounces, pounds, and grains
- Automatic calibration
- Parts counting and percent difference
- Variable integration time
- Stainless steel weighing pan
- Overload protection
- Below-pan weighing
- UL approved
- Spill cover
- DI's 5-year, no-nonsense warranty



SPECIFICATIONS


DI MODEL	DENVER INSTRUMENT CAT NO	CAPACITY	SENSITIVITY	PAN SIZE	REPEATABILITY	LINEARITY	LIST PRICE
XL-410	8515.1	410g	0.001g	4"	0.001g	±0.002g	\$1,395
XL-610	8516.1	610g	0.01g	5.5"	0.01g	±0.01g	\$1,050
XL-1810	8517.1	1810g	0.01g	5.5"	0.01g	±0.01g	\$1,245
XL-3100	8518.1	3100g	0.01g	5.5"	0.01g	±0.01g	\$1,350
XL-6100	8519.1	6100g	0.1g	8.5 x 7"	0.1g	±0.1g	\$1,125
XL-410D	8520.1	410/61g	0.01/0.001g	4"	0.01/0.001g	±0.01/0.002g	\$1,325
XL-3100D	8521.1	3100/450g	0.1/0.01g	8.5 x 7"	0.1/0.01g	±0.1/0.01g	\$1,295

COMMON SPECIFICATIONS

Stabilization time	2 seconds
Net weight	7 lbs (XL-6100/XL-3100 weigh 8 lbs)
Temperature specifications	59-104F (15-40C)
Electrical requirements	115 VAC to 12 VDC transformer included. 230 VAC to 12 VDC transformer available upon request
Dimensions (L x W x H)	10 x 8.5 x 2.5" (25 x 22 x 6 cm)

Denver Instrument Company

XL-SERIES TOPLOADERS



FEATURES

Die-Cast Aluminum
Multiple Weighing Modes
Parts Counting
Stainless Steel Pan
Three Button Operation
Automatic Calibration
Solid-State Electronics
RS-232 Interface
Vacuum-Fluorescent Display
Full-Range Taring
External Power Supply (12V)
In-Use Cover (Standard)
Below-Pan Weighing
Lock Down Device (Optional)
Overload Protection
Three Year Warranty
Made in the U.S.A. (Since 1880)
DEMO EXPRESS
Competition

BENEFITS

Durable, Long Lasting
Versatility
Inventory Control
Non-Corrosive
Easy to Use
Ensures Accurate Readings, Preset
Fast, Reliable Response Time
Transmit Data to Printer/Computer
Bright, Easy to Read
Ability to Subtract Container Weight
Safe, Portable
Spill Protection
Specific Gravity Capability
Theft Protection
Protection Against Abuse
Worry-Free Operation
Over 100 Years of Experience
No Risk Test Drive
METTLER: PB Series
SARTORIUS: BP Series
OHAUS: TP Series
A&D: HF, FX & FY Series

81:

210