


## XL-SERIES TOPLOADING BALANCES

### FEATURES

- RS-232C bi-directional interface
- Multiple weighing units - grams, ounces, pounds, and grains
- Automatic calibration
- Parts counting and percent difference
- Variable integration time
- Stainless steel weighing pan
- Overload protection
- Below-pan weighing
- UL approved
- Spill cover
- DI's 5-year, no-nonsense warranty



### SPECIFICATIONS


DI MODEL	DENVER INSTRUMENT CAT NO	CAPACITY	SENSITIVITY	PAN SIZE	REPEATABILITY	LINEARITY	LIST PRICE
XL-410	8515.1	410g	0.001g	4"	0.001g	±0.002g	\$1,395
XL-610	8516.1	610g	0.01g	5.5"	0.01g	±0.01g	\$1,050
XL-1810	8517.1	1810g	0.01g	5.5"	0.01g	±0.01g	\$1,245
XL-3100	8518.1	3100g	0.01g	5.5"	0.01g	±0.01g	\$1,350
XL-6100	8519.1	6100g	0.1g	8.5 x 7"	0.1g	±0.1g	\$1,125
XL-410D	8520.1	410/61g	0.01/0.001g	4"	0.01/0.001g	±0.01/0.002g	\$1,325
XL-3100D	8521.1	3100/450g	0.1/0.01g	8.5 x 7"	0.1/0.01g	±0.1/0.01g	\$1,295

### COMMON SPECIFICATIONS

Stabilization time	2 seconds
Net weight	7 lbs (XL-6100/XL-3100 weigh 8 lbs)
Temperature specifications	59-104F (15-40C)
Electrical requirements	115 VAC to 12 VDC transformer included. 230 VAC to 12 VDC transformer available upon request
Dimensions (L x W x H)	10 x 8.5 x 2.5" (25 x 22 x 6 cm)

# Denver Instrument Company

## **XL-SERIES TOPLOADERS**



### **FEATURES**

Die-Cast Aluminum  
Multiple Weighing Modes  
Parts Counting  
Stainless Steel Pan  
Three Button Operation  
Automatic Calibration  
Solid-State Electronics  
RS-232 Interface  
Vacuum-Fluorescent Display  
Full-Range Taring  
External Power Supply (12V)  
In-Use Cover (Standard)  
Below-Pan Weighing  
Lock Down Device (Optional)  
Overload Protection  
Three Year Warranty  
Made in the U.S.A. (Since 1880)  
**DEMO EXPRESS**  
Competition

### **BENEFITS**

Durable, Long Lasting  
Versatility  
Inventory Control  
Non-Corrosive  
Easy to Use  
Ensures Accurate Readings, Preset  
Fast, Reliable Response Time  
Transmit Data to Printer/Computer  
Bright, Easy to Read  
Ability to Subtract Container Weight  
Safe, Portable  
Spill Protection  
Specific Gravity Capability  
Theft Protection  
Protection Against Abuse  
Worry-Free Operation  
Over 100 Years of Experience  
No Risk Test Drive  
**METTLER: PB Series**  
**SARTORIUS: BP Series**  
**OHAUS: TP Series**  
**A&D: HF, FX & FY Series**

81:

210