

## DLT Series Portable Balances



# DLT Series Portable Toploading Balances

## The ultimate low-cost balance!

The Denver Light series provides a low-cost alternative for the most basic weighing applications.

**Battery operated.** Powered by a single 9V battery the DLT Series provides the ultimate portability. Purchase the optional transformer for power when near a wall outlet.

**Light-weight.** The DLT Series weighs only 0.9 lbs (0.4 kg). The light-weight but sturdy design is perfect for on-the-go measurements.

**Energy-saving mode.** When powered by battery, the balance will automatically turn itself off after 6 minutes of inactivity.

**Large weighing pan.** The 5.3" x 5.5" pan provides ample space for most loads. The pan can also be removed for quick clean-up.

**Simple three-key operation.** Press the ON/OFF button to cycle power. Dedicated ZERO key tares the balance. The MODE/CAL button toggles between weigh modes and, when held down, initiates the calibration routine.

**Multiple weigh units.** Choose from grams, ounces, pound, pennyweight or troy ounces.



Model	DLT-211	DLT-411	DLT-2100	DLT-4100
Capacity	210g	410g	2100g	4100g
Readability	0.1g	0.1g	1g	1g
Taring Range	0 to 210g	0 to 410g	0 to 2100g	0 to 4100g
Repeatability	0.1g	0.2g	1g	2g
Linearity	0.2g	0.2g	2g	2g
Stabilization Time	3 seconds	3 seconds	3 seconds	3 seconds
Pan Size	5.3 x 5.5" (135 x 140mm)	5.3 x 5.5" (135 x 140mm)	5.3 x 5.5" (135 x 140mm)	5.3 x 5.5" (135 x 140mm)

**Common Specifications:**

Dimensions (L x W x H): 8.5 x 6.3 x 2.0" (217 x 159 x 53mm)  
Net Weight: 0.9 lbs (0.4 kg)  
Operating Temperature: 10°C – 30°C (50°F – 86°F)  
Electrical Requirements: 9V alkaline battery or AC adapter (sold separately)

 **DENVER  
INSTRUMENT**  
Precision laboratory instruments since 1880

6542 Fig St.  
Arvada, Colorado 80004  
800.321.1135  
303.431.7255  
303.423.4831 FAX  
[www.denverinstrumentUSA.com](http://www.denverinstrumentUSA.com)